Carcinoembryonic Antigen

(CEA) IMR Reagent

REF MF-CEA-0061

For *In Vitro* Diagnostic & Professional Use

Intended Use

Use Carcinoembryonic Antigen (CEA) IMR Reagent only with the XacPro-E System (MagQu Co., Ltd.). The reagent is used to quantitatively measure CEA protein in human fluid specimen, such as serum.

Summary & Explanation

Human carcinoembryonic antigen (CEA) (MW = 175,000 – 200,000 g/mol) is a glycoprotein involved in cell adhesion. It is normally produced during fetal development, but the production of CEA stops before birth. Therefore, it is not usually present in the blood of healthy adults, although levels are raised in heavy smokers. Both benign and malignant (harmless and cancerous) conditions can increase the CEA level in blood. The most frequent cancer which causes an increased CEA concentration is cancer of the colon and rectum (Colorectal Cancer, CRC). The normal range for CEA concentration in an adult non-smoker is < 2.5 ng/ml and for a smoker < 5.0 ng/ml in blood.^{1,2}

CEA (ng/mL)	ColoRectal Cancer
> 10	High risk
2.5 (5.0) -10	Low risk
< 2.5 (5.0)	Normal (Smoker)

Principles of Test

The Carcinoembryonic Antigen (CEA) IMR Reagent is designed for rapid quantifying CEA by Immuno Magnetic Reduction (IMR). We conjugated the antibody on the surface of around 50 nm-in-diameter Fe₃O₄ magnetic particles. When the antibodies on the surface bind with CEA, the magnetic particles form clusters. Therefore, the ac susceptibility (X_{ac}) of magnetic particles would be reduced in the adding ac magnetic field. By measuring the reduction of X_{ac} , we can quantify CEA protein in the sample easily, rapidly, and accurately.^{3,4}

Reagents

Carcinoembryonic Antigen (CEA) IMR reagent4 x 1 mL (100 tests)

Storage Conditions & Stability

Storage reagent at 2 ~8 °C, the shelf life is 12 weeks.

Statement of Warnings

1. Do not be frozen.
2. Please keep away from events with strong magnetism.
3. For *in vitro* diagnostic use only.
4. For professional use only.
5. Do not use the reagent when it has taken out form 2~8 °C over 24 hours. Do not use the reagent when it has something precipitated.

Reagent Preparation

Please use the Carcinoembryonic Antigen (CEA) IMR reagents at room temperature (15-30°C).

Specimen Collection & Preparation

1. Collect all blood samples by wearing protective equipment and following universal precautions for venipuncture. Whole blood is collected in an anticoagulant-free serum separation tube (red top).
2. Ensure that complete clot formation has taken place prior to centrifugation.
3. Centrifuge the serum separation tube for 10 minutes at 1,500-2,000 x g to separate the serum from the red blood cells.
4. After centrifugation, serum is taken for measurement immediately or deep frozen (-20°C or colder) until needed. Avoid repeated freezing or thawing.

Procedure

1. Add 60μ L sample to glass testing tube.
2. Add 40μ L magnetic reagent to glass tube, and mix well.
3. Put the glass tube on the measuring site of IMR analyzer.
4. Wait for the end of reaction and read the result.
5. Use the table 1 for converting to the concentration of CEA antigen as reference.

Results

By using XacPro-E, we can get two signals: one is the AC signal before the reaction (X_{ac_0}) and the other is the AC signal after reaction (X_{ac}). Then we can have the IMR (%) through two signal by following function :

$$IMR(\%) = \frac{X_{ac_0} - X_{ac}}{X_{ac_0}}$$

IMR (%), as functions of CEA antigen concentration ϕ_{CEA} are explored and are found to follow the logistic function:

$$IMR(\%) = \frac{A - B}{1 + (\frac{\phi_{CEA}}{\phi_o})^\gamma} + B$$

where A, B, ϕ_o , and γ are fitting parameters. For CEA, A = 1.01, B = 4.25, ϕ_o = 33.35, and γ = 0.76. The concentration of CEA antigen can be available by following equation. And you can find the table for converting to the concentration of CEA as reference.


Fig.1 The IMR standard curve of CEA concentration in serum.

Limitations

1. The analytical range of reagent is from 0.13 to 100 ng/mL. When the specimen with CEA protein > 100 ng/mL is to be determined, carry out the following procedures to obtain the accurate concentration. Dilute the specimen, re-assay, and multiply the assayed CEA protein value by the dilution factor.
2. Reagents should be used before the expiration date printed on the kit label.
3. Data is based upon human serum sample.

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EC REP MedNet GmbH